

43

Surgery  
Chirurgie  
Cirugía  
Chirurgie  
Chirurgia



Couplant Replacement

43-101-15	9 cm / 3½"	Ø 1.5 mm
43-101-20	9 cm / 3½"	Ø 2.0 mm
43-101-25	9 cm / 3½"	Ø 2.5 mm
43-101-35	9 cm / 3½"	Ø 3.5 mm

Tip laterally perforated

43-102-15	9 cm / 3½"	Ø 1.5 mm
43-102-20	9 cm / 3½"	Ø 2.0 mm
43-102-25	9 cm / 3½"	Ø 2.5 mm
43-102-35	9 cm / 3½"	Ø 3.5 mm



Coupland


41-103-16	16.5 cm - 6½"
-----------	---------------



Byrd

43-104-15	21.5 cm / 8½"	Ø 1.5 mm
43-104-25	21.5 cm / 8½"	Ø 2.5 mm
43-104-35	21.5 cm / 8½"	Ø 3.5 mm

Suction cannula with  
finger cut-off and stylet



43-105-15	14 cm / 5½"	Ø 1.5 mm
43-105-20	14 cm / 5½"	Ø 2.0 mm
43-105-25	14 cm / 5½"	Ø 2.5 mm
43-105-30	14 cm / 5½"	Ø 3.0 mm
43-105-40	14 cm / 5½"	Ø 4.0 mm
43-105-50	14 cm / 5½"	Ø 5.0 mm
43-105-60	14 cm / 5½"	Ø 6.0 mm
43-106-15	17 cm / 6¾"	Ø 1.5 mm
43-106-20	17 cm / 6¾"	Ø 2.0 mm
43-106-25	17 cm / 6¾"	Ø 2.5 mm
43-106-30	17 cm / 6¾"	Ø 3.0 mm
43-106-40	17 cm / 6¾"	Ø 4.0 mm
43-106-50	17 cm / 6¾"	Ø 5.0 mm
43-106-60	17 cm / 6¾"	Ø 6.0 mm

Frazier		
43-107-20	12.5 cm / 5"	Ø 2.0 mm
43-107-23	12.5 cm / 5"	Ø 2.3 mm
43-107-27	12.5 cm / 5"	Ø 2.7 mm
43-107-30	12.5 cm / 5"	Ø 3.0 mm
43-107-33	12.5 cm / 5"	Ø 3.3 mm
43-107-37	12.5 cm / 5"	Ø 3.7 mm
43-107-40	12.5 cm / 5"	Ø 4.0 mm
43-107-50	12.5 cm / 5"	Ø 5.0 mm
19 cm / 7½"	30°	
43-108-20	12.5 cm / 5"	Ø 2.0 mm
43-108-23	12.5 cm / 5"	Ø 2.3 mm
43-108-27	12.5 cm / 5"	Ø 2.7 mm
43-108-30	12.5 cm / 5"	Ø 3.0 mm
43-108-33	12.5 cm / 5"	Ø 3.3 mm
43-108-37	12.5 cm / 5"	Ø 3.7 mm
43-108-40	12.5 cm / 5"	Ø 4.0 mm
43-108-50	12.5 cm / 5"	Ø 5.0 mm
19 cm / 7½"	75°	

Fergusson		
43-109-15	11 cm / 4¼"	Ø 1.5 mm
43-109-20	11 cm / 4¼"	Ø 2.0 mm
43-109-25	11 cm / 4¼"	Ø 2.5 mm
43-109-30	11 cm / 4¼"	Ø 3.0 mm
43-109-40	11 cm / 4¼"	Ø 4.0 mm
43-109-50	11 cm / 4¼"	Ø 5.0 mm
18.5 cm / 7¼"		
43-110-15	13 cm / 5"	Ø 1.5 mm
43-110-20	13 cm / 5"	Ø 2.0 mm
43-110-25	13 cm / 5"	Ø 2.5 mm
43-110-30	13 cm / 5"	Ø 3.0 mm
43-110-40	13 cm / 5"	Ø 4.0 mm
43-110-50	13 cm / 5"	Ø 5.0 mm
20.5 cm / 8"		

Suction Tubes / Cleaning brush  
 Saugrohre / Reinigungsbürste  
 Tubos de succión / Cepillo de limpieza  
 Tubes d'aspiration / Brosse de nettoyage  
 Tubi di aspirazione / Pennello di pulizia



Adson

43-111-21 21 cm / 8¼" Ø 4.0 mm



Plester

43-112-10 19.5 cm / 7½" Ø 1.0 mm

43-112-20 19.5 cm / 7½" Ø 2.0 mm

43-112-25 19.5 cm / 7½" Ø 2.5 mm

43-112-30 19.5 cm / 7½" Ø 3.0 mm

For tube with Ø 6 - 9 mm

6.5 mm



43-113-18 18 cm / 7" Ø 5.0 mm

For tube with Ø 6 - 9 mm



Luer


43-114-30 11 cm / 4¼" Ø 3.0 mm

43-114-35 11 cm / 4¼" Ø 3.5 mm

43-114-40 11 cm / 4¼" Ø 4.0 mm

16.5 cm / 6½"

Suction tube, conical



43-115-23

23 cm / 9"



No. 3  
43-116-03  
12 cm / 4<sup>3</sup>/<sub>4</sub>"



No. 3  
43-117-03  
12 cm / 4<sup>3</sup>/<sub>4</sub>"  
cm graduation



No. 3  
43-118-03  
12 cm / 4<sup>3</sup>/<sub>4</sub>"  
cm and inch  
graduation



No. 4  
43-119-04  
13,5 cm / 5<sup>1</sup>/<sub>4</sub>"



No. 4  
43-120-04  
13,5 cm / 5<sup>1</sup>/<sub>4</sub>"  
cm graduation



No. 3R  
43-121-03  
14,5 cm / 5¾"



No. 3R  
43-122-03  
14,5 cm / 5¾"



No. 3R  
43-123-03  
15 cm / 6"



No. 5  
43-124-05  
16 cm / 6¼"



No. 7  
43-125-07  
16 cm / 6¼"



No. 7K  
43-126-07  
12,5 cm / 5"



43-127-15  
15,5 cm / 6"



43-128-13  
13,5 cm / 5¼"




43-129-13  
9,5 cm / 3¾"



Parker  
43-130-14  
14 cm / 5½"




Parker  
43-131-14  
14 cm / 5½"



43-132-06  
6 cm / 2½"  
Scalpel blade remover




43-133-14  
14 cm / 5½"




43-134-13	13 cm / 5"	I	A
43-134-14	14,5 cm / 5¾"	I	A
43-135-13	13 cm / 5"	C	A
43-135-14	14,5 cm / 5¾"	C	A

**XtremeCut**





43-136-14	14,5 cm / 5¾"	I	A
43-137-14	14,5 cm / 5¾"	C	A



43-138-13	13 cm / 5"	I	A
43-138-14	14,5 cm / 5¾"	I	A
43-139-13	13 cm / 5"	C	A
43-139-14	14,5 cm / 5¾"	C	A


**XtremeCut**



43-140-14	14,5 cm / 5¾"	I	A
43-141-14	14,5 cm / 5¾"	C	A





**XtremeCut**




43-142-13	13 cm / 5"	I	△
43-142-14	14,5 cm / 5¾"	I	△
43-143-13	13 cm / 5"	⌋	△
43-143-14	14,5 cm / 5¾"	⌋	△

43-144-14	14,5 cm / 5¾"	I	△
43-145-14	14,5 cm / 5¾"	⌋	△



43-146-14	14,5 cm / 5¾"	I	△	I
43-147-14	14,5 cm / 5¾"	I	△	I
43-148-14	14,5 cm / 5¾"	I	△	I


43-149-14	14,5 cm / 5¾"	I	△	I
43-150-14	14,5 cm / 5¾"	I	△	I
43-151-14	14,5 cm / 5¾"	I	△	I



43-152-12	12 cm / 4¾"	☐	☐
43-153-12	12 cm / 4¾"	☐	☐
43-154-12	12 cm / 4¾"	☐	☐
43-155-12	12 cm / 4¾"	☐	☐



43-156-13	13 cm / 5"	☐	☐
43-157-13	13 cm / 5"	☐	☐



**Deaver**

43-158-14	14 cm / 5½"	☐	☐
43-159-14	14 cm / 5½"	☐	☐



**Kelly**

43-160-16	16 cm / 6¼"	☐	☐
43-161-16	16 cm / 6¼"	☐	☐

**XtremeCut**

43-162-16	16 cm / 6¼"	☐	☐
43-163-16	16 cm / 6¼"	☐	☐

**XtremeCut**




Joseph

43-164-14	14 cm / 5½"	I	△
43-165-14	14 cm / 5½"	J	△



Joseph



43-166-14	14 cm / 5½"	I	△
43-167-14	14 cm / 5½"	J	△



Sanvenero

43-168-14	14 cm / 5½"	I	△
43-169-14	14 cm / 5½"	J	△

**XtremeCut**





Metzenbaum

43-170-11	11.5 cm / 4½"	I	△
43-170-14	14.5 cm / 5¾"	I	△
43-170-18	18 cm / 7"	I	△
43-171-11	11.5 cm / 4½"	U	△
43-171-14	14.5 cm / 5¾"	U	△
43-171-18	18 cm / 7"	U	△

Metzenbaum

43-172-11	11.5 cm / 4½"	I	△
43-172-14	14.5 cm / 5¾"	I	△
43-172-18	18 cm / 7"	I	△
43-173-11	11.5 cm / 4½"	U	△
43-173-14	14.5 cm / 5¾"	U	△
43-173-18	18 cm / 7"	U	△



Metzenbaum

43-174-14	14.5 cm / 5¾"	I	△
43-175-14	14.5 cm / 5¾"	I	△
43-176-14	14.5 cm / 5¾"	U	△
43-177-14	14.5 cm / 5¾"	U	△

Lexer - Fino


43-178-16	16.5 cm / 6½"	I	△
43-179-16	16.5 cm / 6½"	U	△

De Bakey

43-180-17	17.5 cm / 7"	U	△
-----------	--------------	---	---


Arteriotomy Scissors

XtremeCut



Metzenbaum - Fino

43-181-14	14.5 cm / 5¾"	I	A
43-182-14	14.5 cm / 5¾"	C	A



Metzenbaum - Fino

43-183-14	14.5 cm / 5¾"	I	A
43-184-14	14.5 cm / 5¾"	C	A



Metzenbaum - Fino

43-185-14	14.5 cm / 5¾"	I	A
-----------	---------------	---	---

XtremeCut

XtremeCut




Metzenbaum - Fino

43-186-14	14.5 cm / 5¾"	C	A
-----------	---------------	---	---



Metzenbaum - Fino



43-187-14	14.5 cm / 5¾"	I	A
-----------	---------------	---	---





Metzenbaum - Fino

43-188-14	14.5 cm / 5¾"	C	A
-----------	---------------	---	---






43-189-15 15 cm / 6"  







43-190-15 15 cm / 6"    
 Nerve dissecting Scissors







Jameson - Werber  
 43-191-13 13 cm / 5"    







Jameson  
 43-192-15 15 cm / 6"    




43-193-11 11 cm / 4½"    
 43-194-11 11 cm / 4½"  






Toennis - Adson


43-195-17 17 cm / 7"    







Toennis - Adson

43-196-17 17 cm / 7"    


XtremeCut







Toennis - Adson

43-197-17 17 cm / 7"    


XtremeCut







Kilner

43-198-12 12 cm / 4<sup>3</sup>/<sub>4</sub>"    
43-198-15 15 cm / 6"  







Kilner

43-199-12 12 cm / 4<sup>3</sup>/<sub>4</sub>"    
43-199-15 15 cm / 6"  



Reynolds

43-200-15 15 cm / 6"    
43-200-18 18 cm / 7"  

XtremeCut



Iris

43-201-09	9 cm / 3½"		
41-201-10	10.5 cm / 4¼"		
41-201-11	11.5 cm / 4½"		
41-201-12	12.5 cm / 5"		



Iris

43-202-09	9 cm / 3½"		
41-202-10	10.5 cm / 4¼"		
41-202-11	11.5 cm / 4½"		
41-202-12	12.5 cm / 5"		



Iris


43-203-11	11.5 cm / 4½"		
-----------	---------------	--	--

XtremeCut



Iris

43-204-11	11.5 cm / 4½"		
-----------	---------------	--	--



Iris

43-205-11	11.5 cm / 4½"		
-----------	---------------	--	--



Iris

43-206-11	11.5 cm / 4½"		
-----------	---------------	--	--






Fig. 1



Fig. 2




Fig. 3

**Iris**


43-207-11	11.5 cm / 4½"	Fig. 1		
43-208-11	11.5 cm / 4½"	Fig. 2		
43-209-11	11.5 cm / 4½"	Fig. 3		

**XtremeCut**

43-210-11	11.5 cm / 4½"	Fig. 1		
43-211-11	11.5 cm / 4½"	Fig. 2		



43-212-11	11.5 cm / 4½"	Fig. 1		
43-213-11	11.5 cm / 4½"	Fig. 2		



**La Grange - S-Shape**

43-214-11	11.5 cm / 4½"		
-----------	---------------	--	--

**XtremeCut**

43-215-11	11.5 cm / 4½"		
-----------	---------------	--	--



43-216-11	11.5 cm / 4½"		
-----------	---------------	--	--



**Quinby**


43-217-12	12.5 cm / 5"		
-----------	--------------	--	--

**XtremeCut**




Dean

43-218-17 17 cm / 6¾"




Dean

43-219-17 17 cm / 6¾"




Dean

43-220-17 17 cm / 6¾"



Lockin

43-221-16 16 cm / 6¼"



Lockin

43-222-16 16 cm / 6¼"



XtremeCut



Goldman - Fox

43-223-13 13 cm / 5"




Goldman - Fox

43-224-13 13 cm / 5"



Goldman - Fox

43-225-13 13 cm / 5"



XtremeCut



Goldman - Fox

43-226-13 13 cm / 5"



Wagner

43-227-12 12 cm / 4 3/4"



Wagner

43-228-12 12 cm / 4 3/4"





Spencer

43-229-09	9 cm / 3½"	
43-229-11	11 cm / 4½"	
43-229-12	12 cm / 4¾"	



Shortbent

43-230-09	9 cm / 3½"	⌋
-----------	------------	---




Northbent

43-231-09	9 cm / 3½"	⌋
-----------	------------	---



Littauer

43-232-14	14 cm / 5½"	
-----------	-------------	--



Buck

43-233-14	14 cm / 5¾"	
43-233-18	18 cm / 7"	



Eiselsberg

43-234-12	12 cm / 4¾"	
-----------	-------------	--



Universal

43-235-12 12 cm / 4<sup>3</sup>/<sub>4</sub>"



**XtremeCut**

Universal

43-236-12 12 cm / 4<sup>3</sup>/<sub>4</sub>"



Beebe

43-237-10 10 cm / 4"



Beebe

43-238-10 10 cm / 4"



Beebe

43-239-10 10 cm / 4"



Beebe



43-240-12 12 cm / 4<sup>3</sup>/<sub>4</sub>"



Beebe

43-241-12 12 cm / 4<sup>3</sup>/<sub>4</sub>"








Round handle  
 Cutting edge 10 mm

43-242-15	15 cm / 6"	□
43-243-15	15 cm / 6"	⌋
43-244-18	18 cm / 7"	□
43-245-18	18 cm / 7"	⌋
43-246-21	21 cm / 8¼"	□
43-247-21	21 cm / 8¼"	⌋
43-248-23	23 cm / 9"	□
43-249-23	23 cm / 9"	⌋

Flat handle  
 Cutting edge 10 mm

43-250-15	15 cm / 6"	□
43-251-15	15 cm / 6"	⌋
43-252-18	18 cm / 7"	□
43-253-18	18 cm / 7"	⌋
43-254-21	21 cm / 8¼"	□
43-255-21	21 cm / 8¼"	⌋
43-256-23	23 cm / 9"	□
43-257-23	23 cm / 9"	⌋



Round handle  
Cutting edge 10 mm

43-258-15	15 cm / 6"	I
43-259-15	15 cm / 6"	C
43-260-18	18 cm / 7"	I
43-261-18	18 cm / 7"	C
43-262-21	21 cm / 8¼"	I
43-263-21	21 cm / 8¼"	C
43-264-23	23 cm / 9"	I
43-265-23	23 cm / 9"	C

Flat handle  
Cutting edge 10 mm

43-266-15	15 cm / 6"	I
43-267-15	15 cm / 6"	C
43-268-18	18 cm / 7"	I
43-269-18	18 cm / 7"	C
43-270-21	21 cm / 8¼"	I
43-271-21	21 cm / 8¼"	C
43-272-23	23 cm / 9"	I
43-273-23	23 cm / 9"	C



Lister

43-274-11	11 cm / 4¼"	
43-274-14	14 cm / 5½"	
43-274-16	16 cm / 6¼"	
43-274-18	18 cm / 7"	
43-274-20	20 cm / 8"	



Universal

	43-275-15	15 cm / 6"	
	43-276-15	15 cm / 6"	
	43-277-15	15 cm / 6"	
	43-278-15	15 cm / 6"	
	43-279-15	15 cm / 6"	
	43-280-15	15 cm / 6"	
	43-281-15	15 cm / 6"	
	43-282-15	15 cm / 6"	
	43-283-15	15 cm / 6"	
	43-284-18	18 cm / 7"	
	43-285-18	18 cm / 7"	
	43-286-18	18 cm / 7"	
	43-287-18	18 cm / 7"	
	43-288-18	18 cm / 7"	
	43-289-18	18 cm / 7"	
	43-290-18	18 cm / 7"	
	43-291-18	18 cm / 7"	
	43-292-18	18 cm / 7"	




acutegrip




Adson

43-293-12 12 cm / 4¾"




Adson

43-294-12 12 cm / 4¾"





Adson

43-295-12 12 cm / 4¾"




Adson Micro


43-296-12 12 cm / 4¾" 




Adson Micro

43-297-12 12 cm / 4¾" 

acutegrip



Adson Micro

43-298-12 12 cm / 4¾" 

acutegrip



43-299-13	13 cm / 5"
43-299-14	14.5 cm / 5¾"

Standard Tip



43-300-13	13 cm / 5"	I
43-300-14	14.5 cm / 5¾"	I

Narrow Tip



43-301-14	14.5 cm / 5¾"	I
-----------	---------------	---



43-302-13	13 cm / 5"
43-302-14	14.5 cm / 5¾"



43-303-13	13 cm / 5"	I
43-303-14	14.5 cm / 5¾"	I



43-304-16	16 cm / 6½"	I
43-304-18	18 cm / 7"	I

Cross serrated



Cushing

43-305-17 17 cm / 6<sup>3</sup>/<sub>4</sub>"

1:2



Cushing

43-306-17 17 cm / 6<sup>3</sup>/<sub>4</sub>"

acutegrip

0.5 mm



Cushing

43-307-17 17 cm / 6<sup>3</sup>/<sub>4</sub>"



Semken

43-308-12 12 cm / 5"

43-308-15 15 cm / 6"

1:2



Semken

43-309-12 12 cm / 5"

43-309-15 15 cm / 6"

1:2



Gillies


43-310-15 15 cm / 6"



Taylor  
 43-311-17 17 cm / 6¾"



Taylor  
 43-312-17 17 cm / 6¾"




Brophy  
 43-313-20 20 cm / 8"




Brophy  
 43-314-20 20 cm / 8"



Gerald  
 43-315-18 18 cm / 7" 



Gerald  
 43-316-18 18 cm / 7" 



43-317-12 12 cm / 4<sup>3</sup>/<sub>4</sub>"




43-318-15 15 cm / 6"





43-319-15 15 cm / 6"



43-320-15 15 cm / 6" 



43-321-15 15 cm / 6" 



**De Bakey**

43-322-16	16 cm / 6¼"	Fig. 1
43-322-20	20 cm / 8"	Fig. 1
43-323-16	16 cm / 6¼"	Fig. 2
43-323-20	20 cm / 8"	Fig. 2
43-324-16	16 cm / 6¼"	Fig. 3
43-324-20	20 cm / 8"	Fig. 3
43-325-16	16 cm / 6¼"	Fig. 4

**De Bakey**

43-326-16	16 cm / 6¼"	Fig. 5	
43-326-20	20 cm / 8"	Fig. 5	



Fig. 1




Fig. 2



Fig. 3



Fig. 4




Fig. 5



Fig. 6




Fig. 7




Fig. 8



London - College

43-327-01	15 cm / 6"	Fig. 1		
43-327-03	15 cm / 6"	Fig. 3		
43-327-05	15 cm / 6"	Fig. 5		



London - College

43-328-01	15 cm / 6"	Fig. 1		
43-328-02	15 cm / 6"	Fig. 2		
43-328-03	15 cm / 6"	Fig. 3		
43-328-04	15 cm / 6"	Fig. 4		
43-328-05	15 cm / 6"	Fig. 5		
43-328-06	15 cm / 6"	Fig. 6		
43-328-07	15 cm / 6"	Fig. 7		
43-328-08	15 cm / 6"	Fig. 8		



London - College

43-329-01	15 cm / 6"	Fig. 1		
43-329-02	15 cm / 6"	Fig. 2		
43-329-03	15 cm / 6"	Fig. 3		
43-329-04	15 cm / 6"	Fig. 4		
43-329-05	15 cm / 6"	Fig. 5		
43-329-06	15 cm / 6"	Fig. 6		
43-329-07	15 cm / 6"	Fig. 7		
43-329-08	15 cm / 6"	Fig. 8		




Fig. 1



Fig. 2




Fig. 3



Fig. 4




Fig. 5




Fig. 6



London - College

43-330-01	15 cm / 6"	Fig. 1		
43-330-02	15 cm / 6"	Fig. 2		
43-330-03	15 cm / 6"	Fig. 3		
43-330-04	15 cm / 6"	Fig. 4		
43-330-05	15 cm / 6"	Fig. 5		
43-330-06	15 cm / 6"	Fig. 6		
43-330-07	15 cm / 6"	Fig. 7		
43-330-08	15 cm / 6"	Fig. 8		

Meriam

43-331-01	16 cm / 6¼"	Fig. 1		
43-331-03	16 cm / 6¼"	Fig. 3		
43-331-05	16 cm / 6¼"	Fig. 5		

Meriam

43-332-01	16 cm / 6¼"	Fig. 1	
43-332-02	16 cm / 6¼"	Fig. 2	
43-332-03	16 cm / 6¼"	Fig. 3	
43-332-04	16 cm / 6¼"	Fig. 4	
43-332-05	16 cm / 6¼"	Fig. 5	
43-332-06	16 cm / 6¼"	Fig. 6	






Fig. 1



Fig. 2



Fig. 3



Fig. 4




Fig. 5




Fig. 6



**Meriam**

43-333-01	16 cm / 6¼"	Fig. 1	
43-333-02	16 cm / 6¼"	Fig. 2	
43-333-03	16 cm / 6¼"	Fig. 3	
43-333-04	16 cm / 6¼"	Fig. 4	
43-333-05	16 cm / 6¼"	Fig. 5	
43-333-06	16 cm / 6¼"	Fig. 6	

**Flagg**

43-334-01	16 cm / 6¼"	Fig. 1	
43-334-02	16 cm / 6¼"	Fig. 2	
43-334-03	16 cm / 6¼"	Fig. 3	
43-334-04	16 cm / 6¼"	Fig. 4	

**Flagg**

43-335-01	16 cm / 6¼"	Fig. 1	
43-335-02	16 cm / 6¼"	Fig. 2	
43-335-03	16 cm / 6¼"	Fig. 3	
43-335-04	16 cm / 6¼"	Fig. 4	






Fig. 3



Fig. 4

**Landgraf**

43-336-03	15 cm / 6"	Fig. 3	
43-336-04	15 cm / 6"	Fig. 4	




**Perry**

43-337-12	12.5 cm / 5"
-----------	--------------




**Perry**

43-338-12	12.5 cm / 5"	
-----------	--------------	---



Meraim - USA

43-339-17 17 cm / 6<sup>3</sup>/<sub>4</sub>" Fig. 1 





London - College (Modified)

43-340-15 15 cm / 6" 



Crane - Kaplan

43-341-01 15 cm / 6" Fig. 1 

43-341-02 15 cm / 6" Fig. 2 

Pocket markers




Fig. 1

Fig. 2



Gillies

43-342-15 15 cm / 6"



Tjebbes

43-343-15 15.5 cm / 6 1/4"



Surgimax

43-344-15 15 cm / 6"

acutegrip



De Bakey

43-345-16

16 cm / 6½"

1:2



43-346-14

14.5 cm / 5¼"

acutegrip



General

43-347-16

16 cm / 6½"



43-348-14

14.5 cm / 5¼"

acutegrip



Gerald

43-349-18

18 cm / 7"




Gerald

43-350-18

18 cm / 7"





Peet

43-351-11 11 cm / 4½" 




Stieglitz

43-352-13 13 cm / 5¼" 



Stieglitz

43-353-45 13.5 cm / 5¼" 45° 

43-353-90 13.5 cm / 5¼" 90°




Virtus (Ralk)

43-354-15 15 cm / 6" 




Halsted - Mosquito

43-355-12 12.5 cm / 5" 




Halsted - Mosquito

43-356-12 12.5 cm / 5" 




Halsted - Mosquito

43-357-12 12.5 cm / 5" 




Halsted - Mosquito

43-358-12 12.5 cm / 5" 




Kelly

43-359-14 14 cm / 5½" 




Kelly

43-360-14 14 cm / 5½" 




Kelly - Rankin

43-361-16 16 cm / 6¼" 




Kelly - Rankin

43-362-16 16 cm / 6¼" 




Crile

43-363-14 14 cm / 5½" 




Crile

43-364-14 14 cm / 5½" 






Providence - Hospital

43-365-14 14 cm / 5½" 




Providence - Hospital

43-366-14 14 cm / 5½" 




Pean

43-367-14 14 cm / 5½" 




Providence - Pean

43-368-14 14 cm / 5½" 



Providence - Pean

43-369-14 14 cm / 5½" 



Rochester - Ochsner

43-370-14 14 cm / 5½" 



Rochester - Ochsner

43-371-14 14 cm / 5½" 

Towel Clamps  
 Handtuchklemmen  
 Abrazaderas de toalla  
 Pinces à serviette  
 Morsetti per asciugamani



Schaedel

43-372-09 9 cm / 3½"



Backhaus

43-373-09 9 cm / 3½"


43-373-11 11 cm / 4¼"



43-374-11 11 cm / 4¼"

43-374-13 13 cm / 5"

43-374-14 14 cm / 5½"



Allis - Baby

43-375-14 14 cm / 5½"



Allis

43-376-15 15 cm / 6" 4:5

43-377-15 15 cm / 6" 5:6

acutegrip



Babcock

43-378-16 16 cm / 6¼"



Babcock

43-379-16 16 cm / 6¼"




Gross - Maier

43-380-20 20 cm / 8" I

43-381-20 20 cm / 8" C

Fig. 1 = 26 x 10 mm  
 Fig. 2 = 30 x 10 mm



Farabeuf

43-382-12 12 cm / 5¼"




Fig. 1 = 30 x 13 mm  
 Fig. 2 = 34 x 13 mm




Roux



43-383-00


Fig. 1 - 3




43-383-01	14 cm / 5½"	24 x 20 mm / 30 x 28 mm	Fig. 1
43-383-02	15 cm / 6"	29 x 26 mm / 36 x 36 mm	Fig. 2
43-383-03	17 cm / 6¾"	33 x 29 mm / 38 x 43 mm	Fig. 3
43-383-04	18 cm / 7"	29 x 24 mm / 45 x 28 mm	Fig. 4




Senn

43-384-17 17 cm / 6¾"   
 43-385-17 17 cm / 6¾" 



43-386-16 16 cm / 6¼"   
 43-387-16 16 cm / 6¼"   
 43-388-16 16 cm / 6¼" 




43-389-16	16 cm / 6¼"	
43-390-16	16 cm / 6¼"	
43-391-16	16 cm / 6¼"	



Gillies  
43-392-18 18 cm / 7"




Gillies  
43-393-18 18 cm / 7"  
Alar hooklet









Strandell - Stille  
43-394-19 19 cm / 7½"  
43-395-19 19 cm / 7½"



Volkmann  
43-396-21 21.5 cm / 8½"



Volkmann  
43-397-21 21.5 cm / 8½"



**Williger**

43-398-12	12 cm / 5"	
43-399-12	12 cm / 5"	

**Williger**

43-400-12	12 cm / 5"	
43-401-12	12 cm / 5"	

**Williger**

43-402-12	12 cm / 5"	
43-403-12	12 cm / 5"	




Fig. 1




Fig. 2




Fig. 3





Fig. 4

**Langenbeck**

43-404-01	22 cm / 8 <sup>3</sup> / <sub>4</sub> "	30 x 11 mm	Fig. 1
43-404-02	22 cm / 8 <sup>3</sup> / <sub>4</sub> "	30 x 14 mm	Fig. 2
43-404-03	22 cm / 8 <sup>3</sup> / <sub>4</sub> "	30 x 16 mm	Fig. 3
43-404-04	22 cm / 8 <sup>3</sup> / <sub>4</sub> "	40 x 11 mm	Fig. 4



Kocher - Langenbeck

43-405-01	22 cm / 8¾"	40 x 11 mm	Fig. 1
43-405-02	22 cm / 8¾"	55 x 11 mm	Fig. 2




Green

43-406-22	22.5 cm / 8¾"	22 x 22 mm
-----------	---------------	------------




Kocher - Langenbeck

43-407-21	21 cm / 8¼"	23 x 8 mm
-----------	-------------	-----------



Middeldorff

43-408-01	21.5 cm / 8½"	14 x 17 mm	Fig. 1
43-408-02	22 cm / 8¾"	20 x 22 mm	Fig. 2



Obwegeser


43-409-20	20 x 9 mm / 12 mm
43-409-45	45 x 9 mm / 12 mm
43-409-60	60 x 9 mm / 12 mm



Obwegeser


43-410-16	16 cm / 6¼"
-----------	-------------

Chin retractor



Obwegeser

43-411-40	11 x 40 mm
43-411-55	11 x 55 mm
43-411-80	16 x 80 mm



Schaefer

43-412-23	23 cm / 9"
-----------	------------

Cheek retractor





Ward

43-413-15 15.5 cm / 6"



43-414-19 19 cm / 7½"




Bowdler - Henry

43-415-16 16 cm / 6¼"




Kay (Austin)

43-416-08 80 mm



Alm

43-417-01	7 cm / 2¾"	Fig. 1	
43-417-02	10 cm / 4"	Fig. 2	



Weitlaner

43-418-01	13 cm / 5"	Fig. 1	
43-418-02	13 cm / 5"	Fig. 2	





Fig. 1 Fig. 2

Halsey

43-419-01	13 cm / 5"	Fig. 1
43-419-02	13 cm / 5"	Fig. 2

acutegrip



0.4 mm  
Fig. 1 Fig. 2

Halsey

43-420-01	13 cm / 5"	Fig. 1
43-420-02	13 cm / 5"	Fig. 2



0.4 mm  
Fig. 2

Baumgartner

43-421-12	12.5 cm / 5"
43-421-14	14 cm / 5½"

acutegrip



0.5 mm  
Fig. 2

Baumgartner

43-422-12	12.5 cm / 5"
43-422-14	14 cm / 5½"




Fig. 2

Crile - Wood

43-423-15	15 cm / 6"
-----------	------------

acutegrip



0.4 mm  
Fig. 2

Crile - Wood

43-424-15	15 cm / 6"
-----------	------------

acutegrip

0.4 mm



Crile - Murray

43-425-15 15 cm / 6" I



Crile - Murray

43-426-15 15 cm / 6" I



Crile - Murray

43-427-16 16 cm / 6¼"

43-427-18 18 cm / 7"

acutegrip

0.5 mm



Mayo - Hegar

43-428-16 16 cm / 6¼"

43-428-18 18 cm / 7"

acutegrip

0.4 mm



Mayo - Hegar

43-429-18 18 cm / 7" I

acutegrip

0.4 mm



De Bakey

43-430-18 18 cm / 7"

acutegrip

0.4 mm



Ryder

43-431-15	15 cm / 6"
43-431-18	18 cm / 7"

acutegrip

0.3 mm



Ryder

43-432-15	15.5 cm / 6"
43-432-18	18 cm / 7"

acutegrip

0.3 mm



43-433-14	14 cm / 5½"	
43-433-18	18 cm / 7"	



Toennis

43-434-18	18 cm / 7"
-----------	------------

acutegrip

0.5 mm



Toennis

43-435-18	18 cm / 7"
-----------	------------

acutegrip

0.5 mm



Toennis

43-436-20	20 cm / 8"
-----------	------------

acutegrip

0.5 mm




Mathieu

43-437-14	14 cm / 5½"
43-437-17	17 cm / 6¾"



Mathieu - Standard Tip

43-438-14	14 cm / 5½"
43-438-17	17 cm / 6¾"



Mathieu

43-439-14	14 cm / 5½"	
43-439-17	17 cm / 6¾"	

acutegrip

0.5 mm



Mathieu - Narrow Tip

43-440-14	14 cm / 5½"	
43-440-17	17 cm / 6¾"	



Olsen - Hegar

43-441-14	14 cm / 5½"
43-441-18	18 cm / 7"

acutegrip

0.5 mm




Olsen - Hegar

43-442-14	14 cm / 5½"
43-442-18	18 cm / 7"

acutegrip

0.5 mm




Olsen - Hegar

43-443-17 17 cm / 6<sup>3</sup>/<sub>4</sub>"



Olsen - Hegar

43-444-17 17 cm / 6<sup>3</sup>/<sub>4</sub>"



Stevens

43-445-13 13 cm / 5"

43-445-15 15 cm / 6"

acutegrip

0.5 mm




Gillies

43-446-17 17 cm / 6<sup>1</sup>/<sub>4</sub>"



Gillies

43-447-17 17 cm / 6<sup>1</sup>/<sub>4</sub>"




Arruga

43-448-17 17 cm / 6<sup>1</sup>/<sub>4</sub>"

acutegrip

0.4 mm






Arruga


43-449-17 17 cm / 6<sup>1</sup>/<sub>4</sub>"



acutegrip

acutegrip






43-450-14    
14 cm / 5½"





43-451-14    
14 cm / 5½"



43-452-18    
18 cm / 7"



43-453-18    
18 cm / 7"



acutegrip




Castroviejo

43-454-14




14 cm / 5½"



Castroviejo

43-455-14

14 cm / 5½"



Castroviejo

43-456-14




17 cm / 6¾"

Titan  
line




acutegrip


0.3 mm



43-457-14   
14 cm / 5½"



43-458-15     
15 cm / 6"



Friedman - Mini  
43-459-12  
12 cm / 4¾"



Friedman  
43-460-14  
14 cm / 5½"



Friedman  
43-461-14  
14 cm / 5½"



Friedman  
43-462-15  
15 cm / 6"



Cleveland  
43-463-14  
14 cm / 5½"



Cleveland  
43-464-17  
17 cm / 6¾"



**Cleveland**

**43-465-15**

15 cm / 6"

Bone cutting forceps



**Mead**

**43-466-16**

16 cm / 6¼"



**Mead**

**43-467-16**

16.5 cm / 6½"



**Blumenthal**

**43-468-01** 15 cm / 6" Fig. 1

**43-468-02** 15 cm / 6" Fig. 2

**43-468-03** 15 cm / 6" Fig. 3




Fig. 1  
30°

Fig. 2  
45°

Fig. 3  
90°



**Luer**

**43-469-01** 15 cm / 6" Fig. 1

**43-469-02** 15 cm / 6" Fig. 2




Fig. 1

Fig. 2



**Beyer**

43-470-18


18 cm / 7"



**Boehler**

43-471-01 15 cm / 6" Fig. 1


43-471-02 15 cm / 6" Fig. 2



**Ruskin**

43-472-01 18 cm / 7" Fig. 1

43-472-02 18 cm / 7" Fig. 2



**Luer**

43-473-01 17 cm / 6¾" Fig. 1

43-473-02 17 cm / 6¾" Fig. 2

43-473-03 17 cm / 6¾" Fig. 3



43-474-20  
 20 cm / 8"



Dingmann  
 43-475-19  
 19 cm / 7½"




43-476-13  
 13 cm / 5¼"  
 Long ratchet



Adson  
 43-477-13  
 13 cm / 5¼"  
 With drill guide



Hartmann  
43-478-16  
16.5 cm / 6½"




Weil - Blakesley  
43-479-11  
11 cm / 4½" (shaft)  
140°



Mead  
43-480-17  
17 cm / 6¾"  
Ø 27 mm  
320 g



43-481-00  
Plastic  
replacement  
discs



Williger  
43-482-26 24 cm / 9½" Ø 26 mm 300 g  
43-482-28 24 cm / 9½" Ø 28 mm 340 g  
43-482-30 24 cm / 9½" Ø 30 mm 400 g  
Lead filled

Mallets  
Schlägel  
Mazos  
Mallettes  
Mallets



**43-101-00**  
1 pair replacement  
Teflon inserts



**43-484-16** 16.5 cm / 7½" Ø 20 mm 130 g


**43-483-18** 18.5 cm / 7½" Ø 25 mm 220 g



**43-485-18**  
18 cm / 7"  
Lead filled



**Williger**  
**43-486-24**  
24 cm / 9½"  
Lead filled



**Ombredanne**  
**43-487-25** 25 cm / 9¾" 720 g





43-488-21  
21 cm / 8½"  
Riveting hammer  
wide edge



43-489-25  
25 cm / 10"  
Riveting hammer  
wide edge



43-490-21  
21.5 cm / 8½"  
Horn mallet



Fig. 0



Fig. 1



Fig. 2



Fig. 3



Fig. 5

Volkmann

43-491-00	17 cm / 6¾"	Fig. 0
43-491-01	17 cm / 6¾"	Fig. 1
43-491-02	17 cm / 6¾"	Fig. 2
43-491-03	17 cm / 6¾"	Fig. 3
43-491-05	17 cm / 6¾"	Fig. 5



Fig. 00/0



Fig. 0/1



Fig. 1/2

Volkmann

43-492-00	13.5 cm / 5¼"	Fig. 00/0
43-492-01	13.5 cm / 5¼"	Fig. 0/1
43-492-02	13.5 cm / 5¼"	Fig. 1/2



Fig. 00/0




Fig. 0/1



Fig. 1/2


Martini

43-493-00	13.5 cm / 5¼"	Fig. 00/0
43-493-01	13.5 cm / 5¼"	Fig. 0/1
43-493-02	13.5 cm / 5¼"	Fig. 1/2




Williger

43-494-00	16.5 cm / 6½"	Fig. 00	
43-495-00	16.5 cm / 6½"	Fig. 0	
43-495-01	16.5 cm / 6½"	Fig. 1	
43-495-02	16.5 cm / 6½"	Fig. 2	




Williger

43-496-00	16.5 cm / 6½"	Fig. 00	
43-497-00	16.5 cm / 6½"	Fig. 0	
43-497-01	16.5 cm / 6½"	Fig. 1	
43-497-02	16.5 cm / 6½"	Fig. 2	




Williger

43-498-00	16.5 cm / 6½"	Fig. 00	
43-499-00	16.5 cm / 6½"	Fig. 0	
43-499-01	16.5 cm / 6½"	Fig. 1	
43-499-02	16.5 cm / 6½"	Fig. 2	



Williger

43-500-00	16.5 cm / 6½"	Fig. 00	
43-501-00	16.5 cm / 6½"	Fig. 0	
43-501-01	16.5 cm / 6½"	Fig. 1	
43-501-02	16.5 cm / 6½"	Fig. 2	



Chompret  
 43-502-15  
 15 cm / 6"



Williger  
 43-503-00 16 cm / 6¼" Fig. 00/0  
 43-504-00 16 cm / 6¼" Fig. 0/0  
 43-504-01 16 cm / 6¼" Fig. 1/1



Martini  
 43-505-00 16 cm / 6¼" Fig. 00/0  
 43-506-00 16 cm / 6¼" Fig. 0/0  
 43-506-01 16 cm / 6¼" Fig. 1/1




Fig. 1R/1L



Fig. 2R/2L



Fig. 3R/3L




Hemingway

43-507-01	18 cm / 7"	Fig. 1R/1L
43-507-02	18 cm / 7"	Fig. 2R/2L
43-507-03	18 cm / 7"	Fig. 3R/3L

Hemingway

43-508-01	18 cm / 7"	Fig. 1R/1L
43-508-02	18 cm / 7"	Fig. 2R/2L
43-508-03	18 cm / 7"	Fig. 3R/3L



Lucas

43-509-01 17 cm / 6¾" Fig. 1/1

43-509-02 17 cm / 6¾" Fig. 2/2

Mead

43-510-17

17 cm / 6¾"  
 Fig. 5

Molt

43-511-16

16 cm / 6¼"



Miller - Colburn  
43-512-18  
18 cm / 7"  
Fig. 1



Miller - Colburn  
43-513-18  
18 cm / 7"  
Fig. 2



Miller - Colburn  
43-514-18  
18 cm / 7"  
Fig. 3



Miller - Colburn  
43-515-18  
18 cm / 7"  
Fig. 4



Miller - Colburn  
43-516-18  
18 cm / 7"  
Fig. 5




Miller - Colburn  
43-517-18  
18 cm / 7"  
Fig. 1




Seldin  
43-518-18  
18 cm / 7"  
Fig. 12







Miller - Colburn  
43-519-17  
17 cm / 6<sup>3</sup>/<sub>4</sub>"



Miller - Colburn  
43-520-17  
17 cm / 6<sup>3</sup>/<sub>4</sub>"



Miller - Colburn  
43-521-17  
17 cm / 6<sup>3</sup>/<sub>4</sub>"



Miller - Colburn  
43-522-17  
17 cm / 6<sup>3</sup>/<sub>4</sub>"



Ochsenbein  
43-523-16  
16 cm / 6<sup>1</sup>/<sub>4</sub>"



Ochsenbein  
43-524-16  
16 cm / 6<sup>1</sup>/<sub>4</sub>"



Williger  
43-525-17  
17 cm / 6¾"



Williger  
43-526-16  
16 cm / 6¼"



Freer  
43-527-18  
18 cm / 7¼"



Freer  
43-528-19  
19 cm / 7½"



Seldin  
43-529-19  
19 cm / 7½"  
Fig. 22



Seldin  
43-530-19  
19 cm / 7½"  
Fig. 23



Ohl  
43-531-18  
18 cm / 7"



Goldman - Fox  
43-532-18  
18 cm / 7"



Williger  
43-533-13  
13 cm / 5"



Williger  
43-534-16  
16 cm / 6¼"



Molt  
43-535-19  
19 cm / 7½"  
Fig. 2/4





Molt  
43-536-19  
19 cm / 7½"  
Fig. 9





Wakefiled  
43-537-17  
17 cm / 6<sup>3</sup>/<sub>4</sub>"  
Fig. 1W



Freer  
43-538-19   
16 cm / 6<sup>1</sup>/<sub>4</sub>"




Freer  
43-539-19   
16 cm / 6<sup>1</sup>/<sub>4</sub>"



Partsch

43-540-02	13.5 cm / 5¼"	2 mm
43-540-03	13.5 cm / 5¼"	3 mm
43-540-04	13.5 cm / 5¼"	4 mm
43-540-05	13.5 cm / 5¼"	5 mm




Partsch

43-541-02	17 cm / 6¾"	2 mm
43-541-03	17 cm / 6¾"	3 mm
43-541-04	17 cm / 6¾"	4 mm
43-541-05	17 cm / 6¾"	5 mm



Partsch

43-542-02	13.5 cm / 5¼"	2 mm
43-542-03	13.5 cm / 5¼"	3 mm
43-542-04	13.5 cm / 5¼"	4 mm
43-542-05	13.5 cm / 5¼"	5 mm




Partsch

43-543-02	17 cm / 6¾"	2 mm
43-543-03	17 cm / 6¾"	3 mm
43-543-04	17 cm / 6¾"	4 mm
43-543-05	17 cm / 6¾"	5 mm



Lexer - Mini

43-544-08	18 cm / 7"	8 mm
43-544-10	18 cm / 7"	10 mm



Lucas

43-545-03	17 cm / 6¾"	3 mm
43-545-04	17 cm / 6¾"	4 mm
43-546-03	17 cm / 6¾"	3 mm
43-546-04	17 cm / 6¾"	4 mm



43-547-11	11.5 cm / 4½"	Ø 1 mm
43-548-13	13 cm / 5"	Ø 1 mm
43-549-14	14.5 cm / 5¾"	Ø 1 mm



- Fig. 00 Ø 0.7 mm
- Fig. 0 Ø 0.8 mm
- Fig. 1 Ø 0.9 mm
- Fig. 2 Ø 1.1 mm
- Fig. 3 Ø 1.3 mm
- Fig. 4 Ø 1.4 mm
- Fig. 5 Ø 1.5 mm
- Fig. 6 Ø 1.6 mm
- Fig. 7 Ø 1.7 mm
- Fig. 8 Ø 1.8 mm

Bowman

43-550-00	13 cm / 5"	Fig. 00/0
43-550-01	13 cm / 5"	Fig. 0/1
43-550-02	13 cm / 5"	Fig. 1/2
43-550-03	13 cm / 5"	Fig. 2/3
43-550-04	13 cm / 5"	Fig. 3/4
43-550-05	13 cm / 5"	Fig. 4/5
43-550-06	13 cm / 5"	Fig. 5/6
43-550-07	13 cm / 5"	Fig. 6/7
43-550-08	13 cm / 5"	Fig. 7/8

Nickel Silver



Bowman

43-551-00	13 cm / 5"	Fig. 00/0
43-551-01	13 cm / 5"	Fig. 0/1
43-551-02	13 cm / 5"	Fig. 1/2
43-551-03	13 cm / 5"	Fig. 2/3
43-551-04	13 cm / 5"	Fig. 3/4
43-551-05	13 cm / 5"	Fig. 4/5
43-551-06	13 cm / 5"	Fig. 5/6
43-551-07	13 cm / 5"	Fig. 6/7
43-551-08	13 cm / 5"	Fig. 7/8

Nickel Silver



43-552-16

16 cm / 6¼"  
 Wooden  
 Pack of 100 pieces



Frenzel

43-553-16

16.5 cm / 6½"



Bruenings

43-554-19


19 cm / 7½"



Buchwald

43-555-19

19 cm / 7¼"



Hartmann

43-556-15

15 cm / 6"



Wieder

43-557-14

14 cm / 5½" 36 mm


43-557-15

15 cm / 6" 43 mm






Tobold  
43-558-14  
14.5 cm / 5¾"




Tobold  
43-559-15  
15 cm / 6"





Doyen  
43-560-13 13.5 cm / 5¼"  
43-560-16 16 cm / 6¼"




24 / 30 mm



30 / 34 mm





Oringer

43-561-01	Small
43-561-02	Medium
43-561-02	Large



Simplex

43-562-11
11 cm / 4¼"



Surgimax

43-563-11
11 cm / 4¼"



Cawood - Minnesota

43-564-16
16 cm / 6½"




Minnesota

43-565-16
16 cm / 6½"



Cawood - Minnesota


43-566-14
14 cm / 5½"




Branemark  
43-567-15  
15 cm / 6"  
Universal vestibulum  
retractor



Surgimax  
43-568-15  
15 cm / 6"  
Universal vestibulum  
retractor




Maty  
43-569-16 16 cm / 6¼" Children  
43-569-17 17 cm / 6¾" Adult



Haken  
43-570-19  
19 cm / 7½"



Haken  
43-571-20  
20 cm / 8"



Lip  
43-572-17  
17.5 cm / 7"



43-573-14  
14 cm / 5½"



43-574-16  
16 cm / 6¼"



43-575-14  
14 cm / 5½"



Doyen - Jansen (Molt)

43-576-12 12 cm / 4¾"

43-576-14 14 cm / 5½"



Denhardt

43-577-12 12 cm / 4¾"



Roser - Koenig

43-578-16 16 cm / 6¼"

43-578-19 19 cm / 7½"



Heister

43-579-13 13 cm / 5"




Jennings

43-580-09 9 cm / 3½"

43-580-11 11 cm / 4¼"

43-580-13 13 cm / 5"

43-580-15 15 cm / 6"



Dingmann  
43-581-00

14.5 cm / 5<sup>3</sup>/<sub>4</sub>"  
Complete with frame,  
3 tongue depressors,  
2 lateral blades,  
2 tooth plates

Dingmann  
43-582-00


14.5 cm / 5<sup>3</sup>/<sub>4</sub>"  
Only frame



43-583-13  
13.5 cm / 5<sup>1</sup>/<sub>4</sub>"  
65 x 25 mm



43-584-13  
13.5 cm / 5<sup>1</sup>/<sub>4</sub>"  
70 x 29 mm



43-585-13  
13.5 cm / 5<sup>1</sup>/<sub>4</sub>"  
80 x 32 mm



**Kilner - Doughty**

**43-586-00**

Chest support  
Adjustable height


Tongue Holding Forceps  
Zungenholdingzange  
Pinzas de sujeción de la lengua  
Langue maternelle Holding Forceps  
Forbici di tenuta della lingua



**Collin**

**43-587-16**    16 cm / 6¼"

**43-587-18**    18 cm / 7"




**43-589-00**

Rubber inserts (pair)

**Young**

**43-588-16**

16 cm / 6¼"




Langenbeck  
 43-590-00  
 Ø 17 - 21 mm  
 Finger protector



Adams  
 43-591-00  
 Napkin holder



Luniatschek  
 43-592-00  
 Gauze packer



McKesson  
 43-593-01 XS - For babies  
 43-593-02 S - For children  
 43-593-03 M - For adults  
 43-593-04 L - For adults  
 43-593-05 L - For edentulous

Plastic coated metal frame without chain