

49

Impression
Eindruck
Impresión
Impression
Impressione

Stainless Steel Impression Trays, Made of 18/8 Stainless Steel
 Edelstahl Abformlöffel aus Edelstahl 18/8
 Bandejas de impresión de acero inoxidable, hechas de acero inoxidable 18/8
 Plateaux d'impression en acier inoxydable, en acier inoxydable 18/8
 Vassoi di impressione in acciaio inox, in acciaio inossidabile 18/8

49-101-00

Solid Partial
 Set of 10



49-102-01
 U1



49-102-02
 U2



49-102-03
 U3



49-102-04
 U4



49-102-05
 U5



49-103-01
 L1



49-103-02
 L2




49-103-03
 L3




49-103-04
 L4



49-103-05
 L5



Upper	U1	U2	U3	U4	U5
A	62 mm	60 mm	53 mm	54 mm	52 mm
B	83 mm	75 mm	73 mm	69 mm	68 mm



Lower	L1	L2	L3	L4	L5
A	63 mm	61 mm	57 mm	60 mm	53 mm
B	82 mm	81 mm	78 mm	74 mm	71 mm

49-104-00

Porforated Partial
Set of 10



49-105-01
U1



49-105-02
U2



49-105-03
U3



49-105-04
U4



49-105-05
U5



49-106-01
L1



49-106-02
L2




49-106-03
L3




49-106-04
L4



49-106-05
L5



Upper	U1	U2	U3	U4	U5
A	62 mm	60 mm	53 mm	54 mm	52 mm
B	83 mm	75 mm	73 mm	69 mm	68 mm



Lower	L1	L2	L3	L4	L5
A	63 mm	61 mm	57 mm	60 mm	53 mm
B	82 mm	81 mm	78 mm	74 mm	71 mm

Stainless Steel Impression Trays, Made of 18/8 Stainless Steel
 Edelstahl Abformlöffel aus Edelstahl 18/8
 Bandejas de impresión de acero inoxidable, hechas de acero inoxidable 18/8
 Plateaux d'impression en acier inoxydable, en acier inoxydable 18/8
 Vassoi di impressione in acciaio inox, in acciaio inossidabile 18/8

49-107-00

Solid Full Denture
 Set of 8



49-108-01
 XL



49-108-02
 L



49-108-03
 M



49-108-04
 S



49-109-01
 XL




49-109-02
 L




49-109-03
 M



49-109-04
 S



Upper	XL	L	M	S
A	70 mm	67 mm	64 mm	60 mm
B	79 mm	76 mm	70 mm	73 mm



Lower	XL	L	M	S
A	64 mm	64 mm	60 mm	57 mm
B	86 mm	70 mm	73 mm	70 mm

49-110-00

Perforated Full Denture
Set of 8



49-111-01
XL



49-111-02
L



49-111-03
M



49-111-04
S



49-112-01
XL




49-112-02
L




49-112-03
M



49-112-04
S



Upper	XL	L	M	S
A	70 mm	67 mm	64 mm	60 mm
B	79 mm	76 mm	70 mm	73 mm



Lower	XL	L	M	S
A	64 mm	64 mm	60 mm	57 mm
B	86 mm	70 mm	73 mm	70 mm

Stainless Steel Impression Trays, Made of 18/8 Stainless Steel
 Edelstahl Abformlöffel aus Edelstahl 18/8
 Bandejas de impresión de acero inoxidable, hechas de acero inoxidable 18/8
 Plateaux d'impression en acier inoxydable, en acier inoxydable 18/8
 Vassoi di impressione in acciaio inox, in acciaio inossidabile 18/8

49-113-00

Solid
 Set of 6



49-114-01
 U1



49-114-02
 U2



49-114-03
 U3




49-115-01
 L1




49-115-02
 L2



49-115-03
 L3



Upper	U1	U2	U3
A	67 mm	64 mm	57 mm
B	81 mm	71 mm	70 mm



Lower	L1	L2	L3
A	60 mm	57 mm	54 mm
B	81 mm	79 mm	76 mm


49-116-00

Perforated
Set of 6



49-117-01

U1



49-117-02

U2



49-117-03

U3




49-118-01

L1




49-118-02

L2




49-118-03

L3



Upper	U1	U2	U3
A	67 mm	64 mm	57 mm
B	81 mm	71 mm	70 mm



Lower	L1	L2	L3
A	60 mm	57 mm	54 mm
B	81 mm	79 mm	76 mm

Stainless Steel Impression Trays, Made of 18/8 Stainless Steel
 Edelstahl Abformlöffel aus Edelstahl 18/8
 Bandejas de impresión de acero inoxidable, hechas de acero inoxidable 18/8
 Plateaux d'impression en acier inoxydable, en acier inoxydable 18/8
 Vassoi di impressione in acciaio inox, in acciaio inossidabile 18/8

49-119-00

Solid Full Denture
 Set of 8



49-120-01
 XL



49-120-02
 L



49-120-03
 M



49-120-04
 S



49-121-01
 XL




49-121-02
 L




49-121-03
 M



49-121-04
 S



Upper	XL	L	M	S
A	70 mm	67 mm	64 mm	60 mm
B	79 mm	76 mm	70 mm	73 mm



Lower	XL	L	M	S
A	64 mm	64 mm	60 mm	57 mm
B	86 mm	70 mm	73 mm	70 mm

49-122-00

Perforated Full Denture
Set of 8



49-123-01
XL




49-123-02
L




49-123-03
M




49-123-04
S




49-124-01
XL




49-124-02
L




49-124-03
M



49-124-04
S



Upper	XL	L	M	S
A	70 mm	67 mm	64 mm	60 mm
B	79 mm	76 mm	70 mm	73 mm




Lower	XL	L	M	S
A	64 mm	64 mm	60 mm	57 mm
B	86 mm	70 mm	73 mm	70 mm

Stainless Steel Impression Trays, Made of 18/8 Stainless Steel
 Edelstahl Abformlöffel aus Edelstahl 18/8
 Bandejas de impresión de acero inoxidable, hechas de acero inoxidable 18/8
 Plateaux d'impression en acier inoxydable, en acier inoxydable 18/8
 Vassoi di impressione in acciaio inox, in acciaio inossidabile 18/8

English Pattern

49-125-00

Mesh Type
 Set of 8




Upper


49-126-00	Size 0
49-126-01	Size 1
49-126-02	Size 2
49-126-03	Size 3

Lower

49-127-00	Size 0
49-127-01	Size 1
49-127-02	Size 2
49-127-03	Size 3



Upper	Size	A	B	C
	0	67 mm	50 mm	12 mm
	1	72 mm	51 mm	13 mm
	2	78 mm	55 mm	15 mm
	3	80 mm	58 mm	15 mm




Lower	Size	A	B	C
	0	70 mm	50 mm	14 mm
	1	77 mm	51 mm	15 mm
	2	79 mm	55 mm	15 mm
	3	80 mm	60 mm	17 mm

English Pattern

49-128-00

Punched Type
Set of 8




Upper


49-129-00	Size 0
49-129-01	Size 1
49-129-02	Size 2
49-129-03	Size 3

Lower

49-130-00	Size 0
49-130-01	Size 1
49-130-02	Size 2
49-130-03	Size 3



Upper	Size	A	B	C
	0	67 mm	50 mm	12 mm
	1	72 mm	51 mm	13 mm
	2	78 mm	55 mm	15 mm
	3	80 mm	58 mm	15 mm



Lower	Size	A	B	C
	0	70 mm	50 mm	14 mm
	1	77 mm	51 mm	15 mm
	2	79 mm	55 mm	15 mm
	3	80 mm	60 mm	17 mm

Stainless Steel Impression Trays, Made of 18/8 Stainless Steel
 Edelstahl Abformlöffel aus Edelstahl 18/8
 Bandejas de impresión de acero inoxidable, hechas de acero inoxidable 18/8
 Plateaux d'impression en acier inoxydable, en acier inoxydable 18/8
 Vassoi di impressione in acciaio inox, in acciaio inossidabile 18/8

English Pattern

49-131-00

Solid Type
 Set of 8




Upper


49-132-00	Size 0
49-132-01	Size 1
49-132-02	Size 2
49-132-03	Size 3

Lower

49-133-00	Size 0
49-133-01	Size 1
49-133-02	Size 2
49-133-03	Size 3



Upper	Size	A	B	C
	0	67 mm	50 mm	12 mm
	1	72 mm	51 mm	13 mm
	2	78 mm	55 mm	15 mm
	3	80 mm	58 mm	15 mm



Lower	Size	A	B	C
	0	70 mm	50 mm	14 mm
	1	77 mm	51 mm	15 mm
	2	79 mm	55 mm	15 mm
	3	80 mm	60 mm	17 mm

English Pattern

49-134-00

Solid Full Denture
Set of 8



49-135-00

U0



49-135-01

U1



49-135-02

U2



49-135-03

U3



49-136-00

L0



49-136-01

L1




49-136-02

L2




49-136-03

L3



Upper	U0	U1	U2	U3
A	67 mm	72 mm	78 mm	80 mm
B	50 mm	51 mm	56 mm	58 mm





Lower	L0	L1	L2	L3
A	70 mm	77 mm	79 mm	80 mm
B	50 mm	51 mm	55 mm	60 mm

Stainless Steel Impression Trays, Made of 18/8 Stainless Steel
 Edelstahl Abformlöffel aus Edelstahl 18/8
 Bandejas de impresión de acero inoxidable, hechas de acero inoxidable 18/8
 Plateaux d'impression en acier inoxydable, en acier inoxydable 18/8
 Vassoi di impressione in acciaio inox, in acciaio inossidabile 18/8


English Pattern

49-137-00

Perforated Full Denture
 Set of 8



Upper	U0	U1	U2	U3
A	67 mm	72 mm	78 mm	80 mm
B	50 mm	51 mm	56 mm	58 mm



Lower	L0	L1	L2	L3
A	70 mm	77 mm	79 mm	80 mm
B	50 mm	51 mm	55 mm	60 mm

49-140-00

Solid - Impression trays for crown and bridge work
Set of 2

49-143-00

Perforated - Impression trays for crown and bridge work
Set of 2



49-141-01
U1P




49-142-01
L1P



49-144-01
U1P




49-145-01
L1P



49-146-01 44 mm
49-146-02 55 mm
49-146-03 60 mm
49-146-04 65 mm

Perforated
For upper right/lower left jaw



49-148-01 44 mm
49-148-02 55 mm
49-148-03 60 mm
49-148-04 65 mm

Perforated
For upper right/lower left jaw



49-150-00
Perforated
Adjustable



49-151-00
Solid
Adjustable

49-147-01 44 mm
49-147-02 55 mm
49-147-03 60 mm
49-147-04 65 mm

Solid
For upper right/lower left jaw

49-149-01 44 mm
49-149-02 55 mm
49-149-03 60 mm
49-149-04 65 mm

Solid
For upper right/lower left jaw